

Disposafe Suction Catheter

MATERIALS USED:

- Poly Vinyl Chloride, Master Batch.

INDICATIONS:

- Retained secretions causing patient's distress or physiological derangement.
- Ineffective cough & unable to clear secretions spontaneously.
- To stimulate cough.
- To obtain sputum sample.

CONTRAINDICATIONS :

- Product should not be used in patient with known hypersensitivity to any of the materials used.
- Patients with severe bleeding disorder, Irritable airway, recent nasal, oral or esophageal surgery.
- Patients with severe gag reflex, occluded nasal passage, nasal bleeding.

INSTRUCTIONS FOR USE:

1. Check the packing carefully, if packing is found damaged, torn or pierced, discard the piece. Ensure proper selection of size before use.
2. Peel open the pack & remove the device carefully.
3. Check the integrity and patency of the device.
4. Connect the universal funnel connector of suction catheter to the regulated / controlled vacuum source.
5. Gently insert the Suction Catheter distal end-smooth, soft & coned open tip suitably into the oro-pharyngeal cavity /patient's airway.
6. Clear any visible secretion before inserting the suction catheter deep into patient's nare, mouth or artificial airway.
7. Do not force catheter in when experiencing resistance, reinsert catheter if needed.
8. Stimulate cough if it does not occur naturally.
9. Withdraw catheter slightly (1-2mm) prior to applying suction.
10. Apply suction only when withdrawing catheter in a rotating manner.
11. No need to rotate catheter if using an in-line/closed suction system.

INSTRUCTIONS FOR USE

12. Total suction duration should not exceed 15 seconds.
13. Rinse suction catheter if necessary.

WARNINGS:

- This device is intended for single use only.
- Do not expose to heat or direct sunlight.
- The use of this product is restricted to a qualified doctor or a Paramedic.
- The product should be used according to the instructions for use, read the instructions carefully before use.
- The product should not be used for experiment.
- Used product may have potential to biohazards.
- Handle and dispose off in accordance with accepted medical practice and applicable local, state and country laws and regulations.

DISPOSAFE DISCLAIMS ANY RESPONSIBILITY FOR POSSIBLE CONSEQUENCES RESULTING FROM IMPROPER USE.

- The product should not be re-processed.
- Do not clean or re-sterilize.
- The product and its packaging must be visually inspect before use. Improper transport and handling may cause structural and/or functional damage to device or packaging.
- The product is guaranteed sterile & non-pyrogenic, if the package has not been opened or damaged.
- The product should be used immediately after opening the packaging.
- For single use only, re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or infections.
- If there is any change, in expected performance of the device or in case of any malfunction, the device should be immediately removed & sent back to supplier/manufacturer for analysis.

CAUTION:

- Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.



Reference No.



Lot No./Batch No.



Manufacturer



Date Of Manufacturing



Use By / Expiry Date



Do Not Reuse



Do Not Resterilize



Non-Pyrogenic



Do Not Use If Package Is Damaged and Consult Instructions For Use



Caution



Consult Instructions For Use



Temperature Limit



Sterilized By Ethylene Oxide Single Sterile Barrier System



Medical Device



Authorised Representative In The European Community



Disposafe Health and Life Care Limited

Plot No. 1 & 2, Phase-II, Sector-59, Ballabhgarh,

Faridabad-121004, Haryana, India

E-mail : info@disposafe.com

Web : www.disposafe.com

Customer Care No. +91-129-4170201

Mfg. Lie. No. MFG/MD/2019/000106



CEpartner4U B.V., ESDOORNLAAN 13,

3951 DB MAARN, THE NETHERLANDS.

E-mail: office@cepartner4u.com

