

Disosafe Mucus Extractor INSTRUCTIONS FOR USE

MATERIAL USED :

- PP & PVC

INDICATIONS :

- It is used to suck the sputum and secretion in trachea by directly inserted into the throat or by inserted tracheal catheter for anaesthesia.

COMPLICATIONS:

- When used for uterine evacuation, vacuum aspirations is 98% effective in removing all uterine contents.
- Retained products of conception require a second aspiration procedure.
- Abnormal uterine bleeding.
- Infertility.
- Infection of the urinary tract, kidney, or blood.

INSTRUCTION FOR USE:

- Option 1: connect the control valve to the green end & then to the suction force: connect the "suction catheter" to the white end
- Option 2: connect the "control valve" with the "suction catheter" and then to the green end connect the suction force to the white end.
- Removing Mucus : Please your baby on her back and tilt her head up slightly. Squeeze the bulb portion of the Mucus extractor to move the air. Gently insert the tip into your baby's nostril and release the bulb. The suction that is created pulls mucus from the baby's nose into the mucus extractor. Expel the mucus into a tissue. Wipe the tip of the mucus extractor with the tissue and repeat as needed. Clear the other nostril in the same way.
- Cleaning : connect the mucus extractor before and after use to prevent passing germs into your baby's nose. Fill a small container with warm, soapy water and soak the mucus extractor for several minutes. Squeeze the bulb into the water and suck some water into it. Shake well and squeeze the water out. Repeat with clean, warm water. Hand the mucus extractor, tip down. All dry before storing.

Warning :

- After use, this product may be a potential biohazard.
- Handle and dispose off in accordance with accepted.
- Medical practice and applicable local, state and country laws and regulations.
- The use of this product is restricted to a qualified doctor or a paramedic.
- Read instruction before use.
- All connections should be hand tightened, without over-tightening, to ensure the security and to prevent damage to the device.
- After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and country laws and regulation.

MANUFACTURER DISCLAIMS ANY RESPONSIBILITY FOR POSSIBLE CONSEQUENCES RESULTING FROM IMPROPER USE.

- The product should not be reprocessed.
- Visually inspect and carefully check the product and packaging before use. Improper transport and handling may cause structural and/or functional damage to device or packaging.
- The product is guaranteed sterile & non-pyrogenic if the package has not been opened or damaged.
- Do not clean or re-sterile. For single use only. Discard after use.
- The product should be used immediately after opening the packaging.
- Law restricts this device to sale by or on the order of a physician
- The device must not be re-used after its usage, it can be a potential biohazard.
- Improper use can result in serious or fatal illness or injury.

STORAGE CONDITION :

- Do not expose to heat or direct sunlight.

REF

Reference No.

LOT

Lot No./Batch No.



Manufacturer



Date Of Manufacturing



Use By / Expiry Date



Sterilized By Ethylene Oxide
Single Sterile Barrier System



Do Not Reuse



Do Not Resterilize



Non-Pyrogenic



Do Not Use If Package Is Damaged



Caution



Consult Instructions For Use



Temperature Limit



Authorised Representative
In The European Community



Medical Device



Disosafe Health and Life Care Limited

Plot No. 1 & 2, Phase-II, Sector-59, Ballabgarh,
Faridabad-121004, Haryana, India

E-mail : info@disosafe.com,

Web : www.disosafe.com

Customer Care No. +91-129-4170201



CEpartner4U B.V.,
ESDOORNLAAN 13,
3951 DB MAARN, THE NETHERLANDS.

E-mail: office@cepartner4u.com

AW/IFU_1379/A, REV:02, Dt. 28-03-2023

CE 0123