

#### MATERIAL USED :

- PVC, ABS, PP, Latex, Silicon, Isoprene, Nylon mesh and SS

#### INDICATIONS :

- Infusion of I.V. Solutions.
- Intermittent intravenous drug administration under gravity feed through a needle or catheter inserted in vein.
- To maintain hydration and/or correct dehydration if patient is unable to take sufficient volume of oral fluid.

#### CONTRAINDICATIONS :

- Product should not be used in patients with known Hypersensitivity to any of the materials used.
- Administration of high viscous fluids.
- Device should not be used for administration of blood and blood components.
- Device must not be used other than as indicated in intended use.

#### INSTRUCTION FOR USE :

- Check the package sealing integrity to ensure sterility.
- Remove the Infusion Set using aseptic technique.
- Open protective cover on the spike ( closure piercing device).
- Insert spike into the out let of the solution container.
- Fill the drip chamber up to half by squeezing and releasing. Do not fill completely.
- Make sure that flow regulator is in 'on' position.
- Release clamp on IV set to allow fluid to flow through.
- Prime flow regulator with 10-20 ml.
- Re-clamp proximal (input) tube.
- Connect output tube to IV access.
- Set desired flow rate by rotating graduated marker and also by drop counting in drip chamber of I.V. Set

## INSTRUCTIONS FOR USE

and please note the drop size from the packing of I.V. Set flow rate have a tolerance of  $\pm 10\%$ .

- Release clamp to allow fluid flow.
- Observe for adequate flow.

#### WARNING :

- The use of this product is restricted to a qualified doctor or a Paramedic
- The product should be used according to instruction.
- After use, this product may be a potential biohazard. Handle and dispose off in accordance with accepted medical practice and applicable local, state and country laws and regulations.
- If there is any change in expected performance of the device or in case of any malfunction the device should be immediately removed & sent back to supplier for analysis

#### MANUFACTURER DISCLAIMS ANY RESPONSIBILITY FOR POSSIBLE CONSEQUENCES RESULTING FROM IMPROPER USE.

- The product should not be reprocessed.
- Visually inspect and carefully check the product and packaging before use. Improper transport and handling may cause structural and/or functional damage to device or packaging.
- The product is guaranteed non-toxic, sterile and non-pyrogenic if the packages is not opened or damaged.
- Do not clean or re-sterile. For single use only. Discard after use.
- The product should be used immediately after opening the packaging.
- The device must not be re-used after its usage, it can be a potential biohazard.
- Improper use can result in serious or fatal illness or injury.

#### STORAGE CONDITION :

- Do not expose to heat or direct sunlight.



Reference No.



Lot No./Batch No.



Manufacturer



Date of Manufacturing



Use By / Expiry Date



Sterilised By Ethylene Oxide Single Sterile Barrier system



Do Not Reuse



Do Not Resterilize



Non-Pyrogenic



Do Not Use If Package is Damaged & Consult Instructions For Use



Temperature Limitation



Caution



Consult Instructions For Use



Authorised Representative In The European Community



Medical Device

### Disposafe Health And Lifecare Limited

Plot No.1 & 2, Phase-II, Sector-59, Ballabhgarh, Faridabad-121004, Haryana (INDIA)  
E-mail : info@disposafe.com  
Web : www.disposafe.com  
Customer Care No. 91-129-4170201  
Mfg. Lic. No. : MFG/MD/2019/000106



CEpartner4U BV., ESDOORNLAAN 13, 3951 DB MAARN, THE NETHERLANDS.  
E-mail : office@cepartner4u.com

AW/IFU\_1029,Rev:02, Dt. 13-07-2023

