

# INSTRUCTIONS FOR USE

## Luer Lock / cap

### MATERIAL USED :

- H.D.P.E (High Density Poly-Ethylene)

### INDICATIONS :

- To be used as a protective stopper

### CONTRAINDICATIONS:

- Product should not be used in patients with known hypersensitivity to the material used.

### INSTRUCTIONS FOR USE:

- Check visually whether the package is intact.
- Open packaging and remove Injection Stopper from sterile package.
- Connect as desired to the female port of the infusion device in use

### WARNING:

- To be use this product by a qualified doctor or a paramedic.
- Read instructions before use.
- The product should be used according to the instruction for use.

### Manufacturer disclaims any responsibility for possible consequences resulting from improper use.

- The product should not be reprocessed.
- Visually inspect and carefully check the product and packaging before use.
- Improper transport and handling may cause structural and /or functional damage to device or packaging.
- The product is guaranteed non-toxic, sterile and Pyrogen free if the package is not opened or damaged.
- Do not clean or re-sterilize. For single use only. Discard after use.
- Re-use of this device may change its mechanical or biological features and may cause failure, allergic reactions of infections.



We Touch...Lives!

## Disposafe Health and Life Care Limited

Plot No. 1 & 2, Phase-II, Sector-59, Ballabhgarh,

Faridabad-121004, Haryana (INDIA)

E-mail : [info@disposafe.com](mailto:info@disposafe.com),

Web : [www.disposafe.com](http://www.disposafe.com)

Customer Care No. +91-129-4170201



CEpartner4U BV., ESDOORNLAAN  
13, 3951 DB MAARN, THE NETHERLANDS.  
E-mail : [office@cepartner4u.com](mailto:office@cepartner4u.com)



Medical Device



Reference No.



Lot No./Batch No.



Manufacturer



Date Of Manufacturing



Use By / Expiry Date



Sterilized by Ethylene oxide  
Single Sterile Barrier system



Temperature Limit



Do Not Reuse



Do Not Resterilize



Non-Pyrogenic



Do Not Use If Package Is Damaged  
and consult instruction for use.



Caution



Consult Instructions For Use



Authorised Representative  
In The European Community



Keep away from sunlight



AW/IFU\_056,REV.:04, Dt. 02-05-2023